

Introducing

PIATE

Group 24

Biosignal Analysis

Microsoft how-old.net age estimation from pictures

Project Oxford



Lack of research concerning face and voice

Palpitate



heart rate from audio/video

A research project estimating heart rate from video and audio

What are we?

A product with applications in medicine, sports and television

User uploads video/audio file

Relevant features are extracted

How it works

Heart rate is estimated using machine learning on a database

Video is streamed back with superimposed heart rate values

Project Management

Group Management

Research

Results

Meetings

Scrum

Slack

Product

Deadlines

Group Management

Research:

Product:

Audio



Back-end •

Video



Front-end (

Data Preparation - HeartAV

Data Preparation

44 subjects - talking in front of a camera, with heart rate measured.

Times in video must be the same for each ML method

Spectrograms for audio


Green intensity of the face for video


Fail hard


Real Time

Smooth


Face Tracker









Machine Learning







Results: Subject Dependent Audio Models



Results: Subject Independent Audio Models




Results: Video Models



Results

| model type | r | rmse |
|--|-------|------|
| | | |
| Audio: Convolutional - Recurrent network | -0.04 | 1.38 |
| - subject independent | | |
| Audio: Convolutional - Recurrent network | 0.68 | 0.75 |
| - subject dependent | | |
| Audio: Recurrent network - subject inde- | 0.05 | 1.33 |
| pendent | | |
| Audio: Recurrent network - subject de- | 0.62 | 0.81 |
| pendent | | |
| Audio: Convolutional - Recurrent network | NaN | 1.19 |
| with data augmentation | | |
| Video: Convolutional - Recurrent network | 0.52 | 1.18 |
| Video: Convolutional - Recurrent network | 0.48 | 1.28 |
| with data augmentation | | |

Architecture Superimposing, Streaming, Serving



Maximum heart rate

Medical


Use Case

Intensity

Current min/max heart rate

Require age and sex

Design




Attractive and minimal















Appendix

• Maximum Heart Rate (MHR) for men = $203.7 \div (1 + e^{(0.033 \times (Age - 104.3))})$

Maximum Heart Rate (MHR) for women = 190.2 ÷ (1 + e^{(0.0453 × (Age - 107.5))})

